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
Panel de control de programa LAS AF

A screenshot of the LAS AF software interface. The top navigation bar includes 'Configuration', 'Acquire', 'Process', 'Quantify', and 'Analysis'. The 'Acquire' tab is active. The interface is divided into several sections: 'Open projects' and 'Acquisition' on the left; a central control area with 'Load | Save | Roi' and 'Load/Save single setting: User Settings'; a 'Visible' control panel with three sliders for 488, 532, and 635 nm; a 'Specimen' section with a diagram of the optical path; a 'Configuration' panel with 'Objective: Empty 1x/0.00', 'Beamsplitter: DD 488/635', and 'Fluo Turret: Scan-BF'; and a 'Spectral' section with a color scale from 400 to 700 nm and a 'PMT 1' control set to 'ON' with 'Gain [V]: 100.0' and 'Offset [%]: 0.00'.

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The image shows a software interface for microscope configuration. At the top, there are three tabs: 'Configuration' (highlighted in red), 'Acquire', and 'Process'. Below the 'Configuration' tab, there is a grid of icons for different settings: 'IPS', 'Laser Config', 'USB Panel', 'Objective' (highlighted in red), and 'Hardware'. An arrow points from the 'Objective' icon to a window titled 'Objective Attributes'. This window displays a table of objective parameters. The 'Resolution XY(488nm)' and 'Resolution Z(488nm)' values are highlighted with a red box.

Objective Attributes			
Type :	HGX PL APO	Resolution XY(488nm) :	139.43
Magnification :	63	Resolution Z(488nm) :	235.82
Numerical Aperture :	1.4	Free Working Distance :	100
Immersion :	OIL	Focus Depth :	0
Coverglass :	0,17	Focus Offset :	0
Phase Ring :		IC Prisms :	E;E1-P
Technique :		Cond. Prism DIC :	K10
Order Number :		Order Number :	11506187

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